

Lower Anterior, Full Crown Restorations

Patient name _____

Item Checkoff List

- New study cast of approved provisionals. Be sure to capture full extension of ridges and flanges.
- New study cast of opposing arch. Be sure to capture full extension of palate and ridges.
- Face bow of upper arch.
Please face-bow mount upper cast in office to avoid distortion in shipment.
- Full series of photographs of provisionals **without** retractors. Take pictures with patient sitting upright or standing, **not** in the dental chair, teeth slightly apart, to capture the correct plane of occlusion. **These pictures are important.** We need this information to evaluate the occlusal plane, front-to-back progression, smile line, mid-line, and overall esthetics, before beginning case.
- Master impression of lower arch with custom tray. Be sure to include full extension of alveolar ridge contours.
- Bite of lower preps in **maximum intercuspation.**
- Cast of originals to be used as reference for surface texture, **lustre** and characterization.
- Detailed instructions for shading, including pictures of preps (if all ceramic) and photos with shade guides.
- Written detailed instructions with description of work.