

Upper Anterior, Full Crown Restorations

Patient name _____

Item Checkoff List

- New study cast of approved upper provisionals. Be sure to capture full extension of palate, ridges, and flanges.
- Face-bow of approved provisional restorations – **not** the preps. **Please face-bow mount upper cast in office to avoid distortion in shipment.**
- New cast of opposing arch. **Be sure to capture full extension of ridges and flanges.**
- Full series of photographs of provisionals **without** retractors. Take pictures with patient sitting upright or standing, **not** in the dental chair, to capture the correct plane of occlusion. **These pictures are important.** We need this information to evaluate the occlusal plane, front-to-back progression, smile line, mid-line, and overall esthetics, before beginning case.
- Master impressions of upper arch with custom tray. Be sure to include full extension of palate, and full alveolar ridge contours.
- Bite of upper anterior preps in **maximum intercuspation.**
- Mounted casts of originals to be used as reference for surface texture, **lustre** and characterization.
- Detailed instructions for shading, including pictures of preps (if all ceramic) and photos with shade guides.
- Written detailed instructions with description of work.